

DyLight 594 Conjugated AffiniPure Donkey Anti-Mouse IgG (H+L)

Catalog No. BA1148

Size 100 μ l

Ig Type IgG

Concentration

1 mg/ml

Immunogen

Mouse IgG (whole molecular)

Form

Concentrated, Liquid

Storage

At 4°C for 12 months from the date of receipt. Avoid repeated freezing and thawing. Protected from light.

Product Name

DyLight 594 Conjugated AffiniPure Donkey Anti-Mouse IgG (H+L)

Host Donkey

Purity

The antibody was purified from antisera by immunoaffinity chromatography.

Contents

1 mg/ml of dyLight 594 conjugated specific antibody;
0.01 M PBS, 5mg BSA/ml, 0.01% Thimerosal, 50% glycerol.

Fluorophore Amax=693nm; Emax=618nm

Specificity

This dyLight 594 conjugated antibody is specific for mouse IgG and shows no cross-reactivity with human/ bovine/rabbit IgG.

Application

	Suggested Dilution Range
Immunofluorescence	1:100-1:500
Flow Cytometry	1:100-1:500

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.