

## DyLight 488 Conjugated AffiniPure Donkey Anti-Mouse IgG (H+L)

**Catalog No.** BA1145

**Size** 100µl

**Ig Type** IgG

**Concentration**

1 mg/ml

**Immunogen**

Mouse IgG (whole molecular)

**Form**

Concentrated, Liquid

**Storage**

At 4°C for 12 months from the date of receipt. Avoid repeated freezing and thawing. Protected from light.

**Product Name**

DyLight 488 Conjugated AffiniPure Donkey Anti-Mouse IgG (H+L)

**Host** Donkey

**Purity**

The antibody was purified from antisera by immunoaffinity chromatography.

**Contents**

1 mg/ml of dyLight 488 conjugated specific antibody;  
0.01 M PBS, 5mg BSA/ml, 0.01% Thimerosal, 50% glycerol.

**Fluorophore** Amax=462nm; Emax=576nm

**Specificity**

This dyLight 488 conjugated antibody is specific for mouse IgG shows no cross-reactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Rabbit, and Sheep Serum Proteins.

**Application**

	Suggested Dilution Range
Immunofluorescence	1:100-1:500
Flow Cytometry	1:100-1:500

Optimal working dilutions must be determined by end users.

**Application Notes**

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.