

DyLight 550 Conjugated AffiniPure Donkey Anti-Rabbit IgG (H+L)

Catalog No. BA1144

Size 100 μ l

Ig Type IgG

Concentration

1 mg/ml

Immunogen

Rabbit IgG (whole molecular)

Form

Concentrated, Liquid

Storage

At 4°C for 12 months from the date of receipt. Avoid repeated freezing and thawing. Protected from light.

Product Name

DyLight 550 Conjugated AffiniPure Donkey Anti-Rabbit IgG (H+L)

Host Donkey

Purity

The antibody was purified from antisera by immunoaffinity chromatography.

Contents

1 mg/ml of dyLight 550 conjugated specific antibody;
0.01 M PBS, 5mg BSA/ml, 0.01% Thimerosal, 50% glycerol.

Fluorophore Amax=562nm; Emax=576nm

Specificity

This dyLight 550 conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with human/ bovine/rabbit IgG.

Application

	Suggested Dilution Range
Immunofluorescence	1:100-1:500
Flow Cytometry	1:100-1:500

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.