

## **BOSTER BIOLOGICAL TECHNOLOGY Co.,Ltd.**

3942 B Valley Ave, Pleasanton, CA, 94566

Phone: 888-466-3604 Fax: 925-215-2184 Email: boster@bosterbio.com Web: www.bosterbio.com

# DyLight 594 Conjugated AffiniPure Donkey Anti-Rabbit IgG (H+L)

Catalog No. BA1147 Product Name

DyLight 594 Conjugated AffiniPure Donkey Anti-Rabbit IgG (H+L)

Ig Type IgG Host Donkey

Purity Concentration \_\_\_\_

The antibody was purified from antisera by immunoaffinity chromatography.

1 mg/ml

Immunogen Contents

Rabbit IgG (whole molecular) 1 mg/ml of dyLight 594 conjugated specific antibody;

Form 0.01 M PBS, 5mg BSA/ml, 0.01% Thimerosal, 50% glycerol.

Concentrated, Liquid

Fluorophore Amax=693nm; Emax=618nm

Storage

At 4°C for 12 months from the date of receipt. Avoid repeated freezing

and thawing. Protected from light.

### **Specificity**

This dyLight 594 conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with human/ bovine/rabbit IgG.

### **Application**

	Suggested Dilution Range
Immunofluorescence	1:100-1:500
Flow Cytometry	1:100-1:500

Optimal working dilutions must be determined by end users.

#### **Application Notes**

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.