

Biotin Conjugated AffiniPure Donkey Anti-Mouse IgG(H+L)

Catalog No. BA1007

Size 100µl

Ig Type IgG

Concentration

1 mg/ml

Immunogen

Mouse IgG (whole molecular)

Form

Concentrated, Liquid

Storage

At 4°C for 12 months from the date of receipt. Avoid repeated freezing and thawing.

Product Name

Biotin Conjugated AffiniPure Donkey Anti-Mouse IgG(H+L)

Host Donkey

Purity

The antibody was purified from antisera by immunoaffinity chromatography.

Conjugation Biotin

Contents

0.1 mg of biotin conjugated specific antibody;
0.01 M PBS, 0.01% Thimerosal, 50% glycerol.

Specificity

Biotin conjugated antibody is specific for Mouse IgG shows no cross-reactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian, Hamster, Horse, Human, Rabbit and Sheep Serum Proteins.

Application

	Suggested Dilution Range
Immunohistochemistry	1:100-1:200
ELISA	1:5000-1:20000

Optimal working dilutions must be determined by end users.

Application Notes

This antibody must be used with various kinds of enzymes or fluorescence conjugated avidin (streptavidin) in IHC and ELISA.

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.