

BOSTER BIOLOGICAL TECHNOLOGY Co.,Ltd.

3942 B Valley Ave, Pleasanton, CA, 94566

Phone: 888-466-3604 Fax: 925-215-2184 Email: boster@bosterbio.com Web: www.bosterbio.com

Biotin Conjugated AffiniPure Donkey Anti-Rabbit IgG(H+L)

Catalog No. BA1002 **Product Name**

Biotin Conjugated AffiniPure Donkey Anti-Rabbit IgG(H+L) Size 100µl

Ig Type IgG Host Donkey

Purity Concentration

The antibody was purified from antisera by immunoaffinity chromatography. 1 mg/ml

Conjugation **Immunogen Biotin**

Rabbit IgG (whole molecular)

Form 0.1 mg of biotin conjugated specific antibody; Concentrated, Liquid

Contents

0.01 M PBS, 0.01% Thimerosal, 50% glycerol.

Specificity At 4°C for 12 months from the date

of receipt. Avoid repeated freezing

and thawing.

Storage

Biotin conjugated antibody is specific for Rabbit IgG shows no cross-reactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian, Hamster, Horse, Human, Mouse, Rat, Sheep Serum Proteins.

Application

	Suggested Dilution Range
Immunohistochemistry	1:100-1:200
ELISA	1:5000-1:20000

Optimal working dilutions must be determined by end users.

Application Notes

This antibody must be used with various kinds of enzymes or fluorescence conjugated avidin (streptavidin) in IHC and ELISA.

Please centrifuge the vial before use.

Dilute with neutral PBS or TBS.

Prepare working solution on day of use.