

 A Meridian Life Science® Company	MyTaq™ HS	MyTaq™	MyFi™	RANGER	VELOCITY	ACCUZYME™
Properties						
Template Length	≤5kb	≤5kb	≤10kb	≤25kb	≤10kb	≤5kb
Hot-Start	✓	-	✓	○	-	-
High-Fidelity	-	-	✓	-	✓	✓
High Specificity PCR	✓	-	✓	○	-	-
Pre Mix	✓	✓	✓	✓	-	✓
Direct Loading	✓	✓	-	-	-	-
Applications						
Standard PCR	✓	✓	✓	○	○	○
Long PCR	-	-	✓	✓	-	-
Fast PCR	✓	✓	✓	-	✓	-
Multiplex PCR	✓	○	○	-	-	-
GC-Rich Templates	✓	✓	✓	-	-	-
Genotyping	✓	✓	✓	-	-	-
Low-Copy PCR	✓	-	✓	-	✓	✓
Universal Polymerases	✓	-	✓	-	-	-
Bisulfite PCR	-	-	-	-	-	-
TA Cloning	✓	✓	✓	✓	-	-
Blunt-End Cloning	-	-	-	-	✓	✓
Colony PCR	✓	✓	-	-	-	-

 Best choice
  Recommended
  Suitable